NanoBone[®]

Bibliography

Bone Formation in a New Dimension

Bibliografie

Knochenaufbau in neuer Dimension

Important notice:

The package insert provided by us sets out procedures for the application of **NanoBone®**. The following applications may deviate from the package insert. The physician is responsible for the choice of procedure. The manufacturer is not liable should an unsuitable procedure be chosen.

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Canullo L, Patacchia O, Sisti A, Heinemann F

Implant Restoration 3 months after One Stage Sinus Lift Surgery in Severely Resorbed Maxillae: 2-Year Results of a Multicenter Prospective Clinical Study

Clin Implant Dent Relat Res 2010, in press

Objectives: This multicenter prospective study was aimed to clinically evaluate implant behavior inserted in severely resorbed maxillae and restored 3 months after sinus grafting.

Materials and Methods: In three clinical centers, 67 totally rough wide diameter implants were inserted during 30 consecutive sinus lifts. Computed tomography and panoramic analysis were preoperatively requested for each patient. Sinus grafting was performed using a nano-crystalline hydroxyapatite sole bone filler; no membrane was used to cover the buccal window. Preoperative residual bone height ranged between 1–4 mm (mean value: 2.70 mm, standard deviation [SD]: 0.9 mm). Uncovering procedure was carried out following 3 months of healing; 2 weeks later, a definitive restoration was seated using platform switching concept. To monitor stability changes, resonance frequency analysis was performed and implant stability quotient (ISQ) values were collected at the first surgery (baseline, T0), at the abutment connection (T1), and at 2-year follow-up (T2). To measure bone changes, patients underwent panoramic analysis after 2-year follow-up. The image analysis software calculated the grafted bone height changes at level of implant site comparing pre-operative and follow-up panoramic films; the software compensated for eventual radiographic distortion.

Results: Mean ISQ value was 35.7 (SD: 8.8) at baseline, 66.61 (SD: 4.76) at T1, and 77.9 (SD: 4.7) at T2. Statistically significant differences ($p \geq 0.005$) regarding ISQ mean values were found between T1 and T0, as well as between T1 and T2. After 24 months of functional loading, only two implants were lost (cumulative survival rate: 97%). During the same observation period, the mean value of radiographic vertical height of grafted sinus was 13.75 mm (SD = 1.3 mm), with a mean gain of 11 mm.

Conclusions: Within the limits of this study, despite of preoperative residual bone height ranging 1 to 4 mm and absence of the membrane covering the buccal bone wall, maxillary sinus lift restoration 14 weeks after first surgery seems to be a reliable procedure using totally-rough surfaced implants restored using platform switching concept and nano-structured hydroxyapatite as sole bone filler.

Heinemann F, Mundt T, Biffar R, Gedrange T, Götz W

A 3-year clinical and radiographic study of implants placed simultaneously with maxillary sinus floor augmentations using a new nanocrystalline hydroxyapatite

J Physiol Pharmacol 2010, in press

The aims of this case series was to evaluate the success rate of implants and their restorations, the sinus bone graft resorption, and the marginal bone loss around the implants when nanocrystalline HA embedded in a silica matrix was exclusively used as grafting material. In 13 partially edentulous patients of a private practice having missing teeth in the posterior maxilla and a subantral bone height between 3 and 7 mm, 19 sinus augmentations (100% NanoBone®, Artoss, Rostock, Germany) by the lateral lift technique were performed. The implants (Tiolox/Tiologic Implants, Dentaaurum, Ispringen, Germany) were simultaneously placed. After 6 to 9 months 37 implants were restored with fixed dental prostheses. The radiographic bone heights over time were estimated with linear mixed models. The implant success rate was 100% after three years. The mean rates of the marginal bone loss over the first year were higher (mesial: -0.55, distal: -0.51 mm) than the annual rates thereafter (mesial:-0.09 mm, distal: -0.08 mm). The mean rates of changes in the total bone height were neglectable (< 0.2 mm) and not significant. The prosthodontic and esthetic evaluation revealed a successful outcome. Within the limits of this clinical report it may be concluded that maxillary sinus augmentation using 100% nanocrystalline HA embedded in a silica matrix to support implants is a reliable procedure.

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Modern bone substitute materials (BSM) have to meet numerous structural and biological requirements. Accordingly, morphological in vitro analysis and histomorphometric ex vivo investigations are of great significance to estimate BSM biocompatibility. Aim of the study was a respective evaluation of a modern alloplastic BSM (NanoBone®). Structural in vitro analysis of the native BSM was carried out by electron microscopy and microcomputed tomography (μ -CT) with special regard to porosity. 14 months after maxillary sinus augmentation with NanoBone® and collected autologous bone particles in one individual patient case, a representative trephine biopsy out of the augmentation volume was histomorphometrically analysed employing conventional histology and μ -CT. Volume ratios of newly formed bone and remaining BSM particles were calculated via assessment of 2-D phase distribution of tissue density. In vitro investigation of the BSM showed a chiselled macrostructure with a total porosity of $> 65\%$ as well as a high ratio of pores $> 250\ \mu\text{m}$, which were almost exclusively localized interparticulary. Histomorphometric analysis of the trephine biopsy revealed a good bony integration of the BSM with evidence of BSM resorption and replacement by vital bone tissue after 14 months. The volume ratio of newly formed bone was 37% . The presented methods for pre-clinical and clinical evaluation of modern BSM complement one another in a reasonable manner.

Mertens C, Steveling H

Use of Synthetic Bone Blocks as an Alternative to Autologous Bone Block Grafts

Implants, International Magazine of Oral Implantology 2009; 4:30-32

In modern implantology, correct three-dimensional positioning of implants, as well as sufficient bone material are of great importance in order to reach satisfactory and predictable results. Resorption processes, traumatic tooth losses or chronic inflammatory processes such as chronic periodontal diseases, however, often result in severe reduction of bone material. If affected areas are intended to serve as implant beds, augmentation will often be required during the same or in a previous intervention. While autologous bone is still considered to be the gold standard, bone substitute materials have proven successful particularly in cases of rather small defects. Their use may decrease patient's morbidity, shorten treatment duration and reduce treatment costs. However, if the defect exceeds a certain size, autologous bone grafts will have to be used, usually in the form of blocks. Intraoral bone

removal poses the problem of limited availability. Extraoral donor sites, however, require treatment under general anesthesia or under in-patient conditions, which is why patients frequently reject this type of surgery. In particular in cases of edentulism in the molar and premolar region, patients tend to prefer fixed dental prostheses, however, the problem of a significantly narrowed alveolar ridge often occurs in the molar area of the mandible.

The use of the **NanoBone® | block** (Artoss, Germany) constitutes a possible alternative to autologous bone blocks. The nanocrystalline material, that has already proven reliable in many trials in a particulate form, has been available on the market in the form of blocks for a short time. Preclinical trials using animal models have shown high rates of bone formation within a relatively short period of time. The following follow-up observation was initiated to find out whether the bone substitute material used in the form of blocks proves successful as a possible alternative to autologous bone.

The nanocrystalline blocks used constitute a possible alternative to autologous bone blocks. The block provides a sufficient primary stability to be used safely for augmentation. The clinical procedure, however, differs from the use of e.g. autologous blocks removed from the retromolar space. The special structure of the block provides for the complete osseointegration of the augmentation material and thus for a sufficient gain in volume for safe implantation

Meier J

Experiences with the nanostructured bone substitute NanoBone™ in particular and block form: Prospective histological and clinical trial with 3 years follow-up

EAO 18th International Meeting Monaco 2009, Poster 245

This study presents the results of different augmentation procedures using the new and nanostructured bone substitute (BS) **NanoBone™** with special regard to the histologic features and demonstrates that former therapy protocols can be changed to remarkable shorter healing periods which can be carried out with reliable results. The structural changes were analysed histologically and the cellular ingrowth of bone forming cell lines and blood vessels could be verified. Based on 86 sinus floor elevations (SFE) and 75 lateral augmentations (LAT) performed on average 3 or more years ago there is no measurable difference in bone height and dimension. Histomorphometry of SFE samples showed about 40 vol.% of de novo bone formation after only 2-3 months which must be compared to other BS.

The preliminary results following augmentations with **NanoBone Blocks™** are encouraging and suggest that this might be a way to abandon the transplantation of bone blocks of other origin.

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Early matrix change of a nanostructured bone grafting substitute in the rat.

J Biomed Mater Res B Appl Biomater. 2009 Nov; 91(2):692-9

A nanocrystalline bone substitute embedded in a highly porous silica gel matrix (**NanoBone®**) has previously been shown to bridge bone defects by an organic matrix. As the initial host response on the bone graft substitute might be a determinant for subsequent bone formation, our present purpose was to characterize the early tissue reaction on this biomaterial. After implantation of 80 mg of **NanoBone®** into the adipose neck tissue of a total of 35 rats, grafts were harvested for subsequent analysis at days 3, 6, 9, 12, and 21. The biomaterial was found encapsulated by granulation tissue which partly penetrated the implant at day 3 and completely pervaded the graft at day 12 on implantation. Histology revealed tartrate-resistant acid phosphatase (TRAP)-positive giant cells covering the biomaterial. ED1 (CD68) immunopositivity of these cells further indicated their osteoclast-like phenotype. Scanning electron microscopy revealed organic tissue components within the periphery of the graft already at day 9, whereas the central hematoma region still presented the silica-surface of the biomaterial. Energy dispersive X-ray spectroscopy further demonstrated that the silica gel was degraded faster in the peripheral granulation tissue

than in the central hematoma and was replaced by organic host components by day 12. In conclusion, the silica gel matrix is rapidly replaced by carbohydrate macromolecules. This might represent a key step in the process of graft degradation on its way toward induction of bone formation. The unique composition and structure of this nanoscaled biomaterial seem to support its degradation by host osteoclast-like giant cells.

Reichert C, Al-Nawas B, Smeets R, Kasaj A, Götz W, Klein MO

In vitro proliferation of human osteogenic cells in presence of different commercial bone substitute materials combined with enamel matrix derivatives

Head & Face Medicine 2009, 5:23

Background: Cellular reactions to alloplastic bone substitute materials (BSM) are a subject of interest in basic research. In regenerative dentistry, these bone grafting materials are routinely combined with enamel matrix derivatives (EMD) in order to additionally enhance tissue regeneration.

Materials and methods: The aim of this study was to evaluate the proliferative activity of human osteogenic cells after incubation over a period of seven days with commercial BSM of various origin and chemical composition. Special focus was placed on the potential additional benefit of EMD on cellular proliferation.

Results: Except for PerioGlas®, osteogenic cell proliferation was significantly promoted by the investigated BSM. The application of EMD alone also resulted in significantly increased cellular proliferation. However, a combination of BSM and EMD resulted in only a moderate additional enhancement of osteogenic cell proliferation.

Conclusion: The application of most BSM, as well as the exclusive application of EMD demonstrated a positive impact on the proliferation of human osteogenic cells in vitro. In order to increase the benefit from substrate combination (BSM + EMD), further studies on the interactions between BSM and EMD are needed.

Punke C, Zehlicke T, Boltze C, Pau HW

Investigation of a new highly porous hydroxyapatite matrix for obliterating open mastoid cavities - application in guinea pigs bulla

Laryngorhinotologie 2009 Apr; 88(4):241-6

Background: Many different techniques for obliterating open mastoid cavity have been described. The results after the application of alloplastic materials like Hydroxyapatite and Tricalciumphosphate were poor due to long-lasting resorption. Extrusion of those materials has been described. We investigated the applicability of a new high-porosity ceramic for obliterating large open mastoid cavities and tested it in an animal model (bulla of guinea pig).

Methods: A highly porous matrix (NanoBone®) bone-inductor fabricated in a sol-gel-technique was administered unilaterally into the opened bullae of 30 guinea pigs. In each animal the opposite bulla was filled with Bio-Oss, a bone substitute consisting of a portion of mineral bovine bone. Histological evaluations were performed 1, 2, 3, 4, 5 and 12 weeks after the implantation.

Results: After the initial phase with an inflammatory reaction creating a loose granulation tissue, we observed the formation of trabecular bone within the fourth week in both groups. From the fifth week on we found osteoclasts on the surface of NanoBone® and Bio-Oss with consecutive degradation of both materials.

Conclusion: In our animal model study we found beneficial properties of the used bone-inductors NanoBone® and Bio-Oss for obliterating open mastoid cavities.

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A prospective study—Histological results after 3 months of healing
Implants, International Magazine of Oral Implantology 2009; 2:24-27

Background: The aim of this prospective study was to evaluate tissue composition of augmented maxillary sinus floor 3 months after using of a nano-crystalline hydroxyapatite bone substitute. Histological analysis and bone-to-implant contact (BIC) assessment between the grafting material and inserted miniimplant were achieved.

Methods: Five patients (2 men and 3 women) in need for fixed implant-supported prosthesis in the posterior maxillae were consecutively recruited for the present study. Preoperatively, computerized tomography and digital panoramic examinations were acquired for antral anatomy evaluation. A rectangular or oval-shaped osteotomy was then prepared on the lateral aspect of the alveolar ridge under copious normal saline irrigation. The resulted detached “window” was elevated medially and apically while simultaneously reflecting the sinus membrane. NanoBone® mixed with antibiotic solution was placed incrementally at the superior aspect of the sinus and against the medial aspect of the grafted compartment created in the sinus cavity. A mini screw for osteosynthesis of 1.2 mm diameter and 13 mm in length was then positioned to maintain the space opened. After a 3 month healing period, a bioptical core containing the mini-implant was retrieved using a 3 mm trephine bur. In the same surgical step, implants were inserted. After 3 months of submerged healing, implants were restored.

Results: After 3 months of healing, varying amounts of newly formed bone were found through the specimens. From the histomorphometric analysis, NanoBone® residuals accounted for 47.35 % ± 5.20 % of the extracted bone volume, marrow spaces presented 19.30 % ± 3.20 % and bone occupied 33.35 % ± 4.1 % (new bone: 22.23 % ± 4.10 %, and native bone: 11.12 % ± 4.20 %). Well-mineralized regenerated bone with lamellar parallel-fibred structure and Haversian systems surrounded the residual NanoBone® particles. Mean BIC was 17.75 % ± 2.9 %. No connective tissue was observed at the implant boundary surface.

Conclusion: Within the limits of this clinical prospective study, it can be concluded that nano-crystalline hydroxyapatite bone substitute showed good histological outcomes for augmenting maxillary sinus floor in critical bone volume conditions. Furthermore, the absence of covering membrane and 3-month healing period could clinically demonstrate the potential of this grafting biomaterial. In such a critical condition the use of a rough-surfaced mini-implant showed BIC values supposed to be effective also in case of functional loading.

Stübinger S, Ghanaati S, Orth C, Hilbig U, Saldamli B, Biesterfeld S, Kirkpatrick CJ, Sader RA

Maxillary sinus grafting with a nano-structured biomaterial: preliminary clinical and histological results.

Eur Surg Res. 2009;42(3):143-9. Epub 2009 Jan 29.

Background: In this study the potential of a new and entirely synthetic, nano-structured hydroxyapatite-based biomaterial for sinus floor augmentation is evaluated.

Methods: 20 sinus floor elevations were carried out in a total of 20 patients. After a healing period of 6 months, in 10 cases cylinder-shaped bone biopsies were taken from the augmented maxillary region using trephine burs.

Results: The healing period progressed without any complications. General and specific histological analysis of the bone biopsies showed a high osteoclast activity at the margin of the biomaterial which was well integrated into the newly formed bone.

Conclusion: This study demonstrates that new trabecular bone is formed after grafting with the nanocrystalline bone substitute after 6 months. Ongoing histomorphological studies are necessary to quantify the biomaterial-bone ratio and the exact amount of newly built bone in the augmented cavity after 6 months.

Canullo L, Dellavia C

Sinus lift using a nano-crystalline hydroxyapatite silica gel in severely resorbed maxillae: histological preliminary study

Clin Implant Dent Relat Res 2009 Oct; 11 Suppl 1:e7-13. Epub 2009 Feb 13

Objectives: The aim of this prospective study was to evaluate histologically, radiographically and clinically a new nano-structured hydroxyapatite in maxillary sinus floor grafting in severely resorbed maxillae. **Materials and methods:** A total of 33 totally length micro-textured implants were placed during 16 consecutive sinus lift. No membrane was used to close the buccal window. Preoperative residual bone level ranged between 1-3 mm (mean value of 2.03 mm). After 3-4 months of healing, definitive restorations were seated using platform switching concept. **Results:** After 18 months of functional loading, no implant was lost. During the same observation period, the mean value of radiographic vertical height of grafted sinus floor was 14.2 mm (SD=0.505 mm) and the mean value of grafted bone filler absorption rate was 4.4 % (SD=0.456 mm). Histological analysis showed significant new bone formation, and remodeling of the grafted material. In the cores obtained at 5 months, regenerated bone, residual NanoBone® and bone marrow occupied respectively 48±4.63 %, 28±5.33 % and 24±7.23 % of the grafted volume. In the specimens taken 3 months after grafting, mean new bone was 8±3.34 %, mean NanoBone® was 45±5.10 % and mean bone marrow was 47±6.81 % of the bioptical volume. **Conclusions:** Within the limits of this study, it was concluded that grafting of maxillary sinus using nano-structured hydroxyapatite as only bone filler is a reliable procedure also in critical conditions.

Abshagen K, Schrodi I, Gerber T, Vollmar B

In vivo analysis of biocompatibility and vascularization of the synthetic bone grafting substitute NanoBone®

J Biomed Mater Res A 2009 Nov; 91(2):557-66

One of the major challenges in application of bone substitutes are adequate vascularization and biocompatibility of the implant. Thus, the temporal course of neovascularization and the microvascular inflammatory response of implants of NanoBone® (fully synthetic nanocrystalline bone grafting material) were studied in vivo by using the mouse dorsal skinfold chamber model. Angiogenesis, microhemodynamics and leukocyte-endothelial cell interaction were analyzed repetitively after implantation in the center and in the border zone of the implant up to 15 days. Both NanoBone® granules and plates exhibited high biocompatibility comparable to that of cancellous bone, as indicated by a lack of venular leukocyte activation after implantation. In both synthetic NanoBone® groups, signs of angiogenesis could be observed even at day 5 after implantation whereas granules showed higher functional vessel density compared with NanoBone® plates. The angiogenic response of the cancellous bone was markedly accelerated in the center of the implant tissue. Histologically, implant tissue showed an ingrowth of vascularized fibrous tissue into the material combined with an increased number of foreign-body giant cells. In conclusion, NanoBone®, particularly in granular form, showed high biocompatibility and high angiogenic response, thus improving the healing of bone defects. Our results underline, that beside the composition and nanostructure, also the macro-structure is of importance for the incorporation of the biomaterial by the host tissue.

Kruse A, Jung RE, Nicholls F, Zwahlen RA, Hämmerle CHF, Grätz KW, Weber FE

Comparison of synthetic HA/SiO₂ matrix and bovine derived HA

EAO 17th International Meeting Warsaw 2008, Poster 243

Background: The substitution of autologous bone with synthetic materials for the treatment of bone defects is still a challenge. Calcium phosphate salts, like hydroxyapatite, are often used to develop synthetic bone substitutes since they are main constituent of natural bone material. During synthesis, most synthetic bone substitutes are sintered

yielding in a more compact and less porous material, where osteoconductivity can be reduced. **NanoBone®**, is a non-sintered nanocrystalline hydroxyapatite embedded in a high porous silica gel matrix. In order to guarantee a high osteoinductive property and a biodegradability, the granula are loosely packed and present a porosity >50 %.

Purpose: The aim of the present study was to test whether or not a synthetic hydroxyapatite/silica oxide based bone substitute material (**NanoBone®**) enhances bone regeneration compared to a xenogenic hydroxyapatite based bone substitute material (**BioOss®**) or empty control sites.

Material and methods: A rabbit calvarial defect model was used to compare the different bone substitute materials. The handling characteristics of both materials was similar. The samples were embedded Goldner Trichrome stained and middle sections were used for evaluation. In none of the sections any signs of inflammation was detectable.

Results: Bone tissue in the defect: The results of the histomorphometric analysis revealed a significant difference between the percentages between bone formed in the empty and the synthetic hydroxyapatite/silica oxide based granules group. **Bone bridging:** Bony bridging is the percentage of the defect where new bone has occurred. The box-plot shows median, the standard deviation and the 95 %-confidence interval (red box). The P values determined by an ANOVA according the Fisher least significant difference Post hoc procedure showed a highly significant increase in bone bridging when the untreated defect group was compared to the group treated with synthetic hydroxyapatite/silica oxide based granules ($P=0.032$). When these two groups were compared by a paired t-test the difference was still significant ($P=0.067$). Both materials show an excellent bone integration.

Conclusion: Compared to empty defects: significantly more bone forms if the defects are treated with synthetic hydroxyapatite/silica oxide granules (**NanoBone®**) significantly more of the defect is bridged by bone when synthetic hydroxyapatite/silica oxide granules are applied (**NanoBone®**). In this in vivo model system no significant difference is seen between hydroxyapatite/silica oxide granules (**NanoBone®**) and xenogenic hydroxyapatite based material (**BioOss®**).

Lenz S, Kirchoff M, Gerber T

Enhanced osseointegration of implants with a nanostructured bioactive coating

EAO 17th International Meeting Warsaw 2008, Poster 391

Objectives: In this study we tried to investigate whether it is possible to use the properties of the bone grafting material **NanoBone®** for coating of dental implants to improve their osseointegration.

Material and methods: The implants (group A: Semados®, sand blasted surface, group B ixx2®, sand blasted and acid etched surface) were coated with a silica matrix covering nanocrystalline hydroxyapatite by sol-gel technique. The implants showed differences in screw thread and roughness. Coated (n=18) and uncoated (n=18) implants were inserted in the frontal bone of 8 minipigs. Specimens were excised after 2, 4 and 6 weeks and processed according to the sawing and grinding technique. The bone to implant contact (BIC) was measured by semiautomatic software.

Results: All coated implants showed a higher rate of BIC compared to the uncoated implants. The mean percentage of BIC for coated implants group A was 60.2 %-2 weeks, 66.6 %-4 weeks, and 74.5 %-6 weeks. The uncoated implants of this group reached 57.0 %-2 weeks, 61.3 %-4 weeks and 64.4 %-6 weeks. In group B the BIC was 73.4 %-2 weeks, 70.6 %-4 weeks and 78.0 % for the coated ones. The uncoated implants in this group reached a BIC of 68.5 %-2 weeks, 60.9 %-4 weeks and 45.8 %-6 weeks.

Conclusion: The applied coating of implants enhances the BIC. Earlier loading of such modified implants can be considered.

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This prospective study evaluated the structural changes (modeling and remodeling) as well as the biodegradation of the new bone grafting material NanoBone® based on clinical and histologic investigation. Sinus floor elevations were performed on 17 patients using a two-stage protocol when the subantral bone height was less than 5 mm. The 43 bone samples were collected during implant placement, which was carried out after healing periods of 8–11 weeks (group I) or 12–15 weeks (group II), and subjected to undecalcified tissue processing by applying a hard specimen cutting-grinding technique. The clinical findings showed a solid ossification with bone qualities of D1 or D2, that could be verified in the histologic sections showing impressive hyperostosis. The resorption of NanoBone® and the de novo bone formation took place simultaneously, similar to the processes following transplantation of autogenous cancellous bone. Compared with other bone substitutes, we observed an accelerated organization and new bone formation that, after only 3 months, yielded a solid bony layer for primary stable implant placement in the augmented maxillary sinus. Early implantation and functional loading stimulated the new bone and prevent a loss.

Götz W, Gerber T, Michel B, Lossdörfer S, Henkel KO, Heinemann F

Immunohistochemical characterization of nanocrystalline hydroxyapatite silica gel (NanoBone®) osteogenesis: A study on biopsies from human jaws

Clin Oral Impl Res 2008; 19;1016-26

Bone substitute biomaterials may be osteogenic, osteoconductive or osteoinductive. To test for these probable characteristics in a new nanoporous grafting material consisting of nanocrystalline hydroxyapatite embedded in a porous silica gel matrix (NanoBone®), applied in humans, we studied biopsies from 12 patients prior to dental implantation following various orofacial augmentation techniques with healing times of between 3.5 and 12 months. Sections from decalcified specimens were investigated using histology, histochemistry (PAS, alcian blue staining, TRAP) and immunohistochemistry with markers for osteogenesis, bone remodelling, resorption and vessel walls (alkaline phosphatase, bone morphogenetic protein-2, collagen type I, ED1, osteocalcin, osteopontin, runx2, vWF). Histologically, four specific stages of graft transformation into lamellar bone could be characterized. During early stages of healing, bone matrix proteins were absorbed by Nanobone® granules, forming a proteinaceous matrix, which was invaded by small vessels and cells. We assume that the deposition of these molecules promotes early osteogenesis in and around NanoBone® and supports the concomitant degradation probably by

osteoclast-like cells. TRAP-positive osteoclast-like cells were localized directly on the granular surfaces. Runx2-immunoreactive pre-osteoblasts, which are probably involved in direct osteogenesis forming woven bone which is later transformed into lamellar bone, were attracted. Graft resorption and bone apposition around the graft granules appear concomitantly. We postulate that NanoBone® has osteoconductive and biomimetic properties and is integrated into the host's physiological bone turnover at a very early stage.

Punke C, Zehlicke T, Boltze C, Pau H-W

Experimental Studies on a New Highly Porous Hydroxylapatite Matrix for Obliterating Open Mastoid Cavities

Otol Neurotol 2008 Sep; 29(6):807-11

Objective: In an initial preliminary study, the applicability of a new high-porosity hydroxyapatite (HA) ceramic for obliterating large open mastoid cavities was proven and tested in an animal model (bulla of guinea pig).

Study Design: Experimental study.

Methods: NanoBone®, a highly porous matrix consisting of 76 % hydroxylapatite and 24 % silicone dioxide fabricated in a sol-gel technique, was administered unilaterally into the opened bullae of 30 guinea pigs. In each animal, the opposite bulla was filled with Bio-Oss, a bone substitute consisting of a portion of mineral bovine bone. Histologic evaluations were performed 1, 2, 3, 4, 5 and 12 weeks after the implantation.

Results: After an initial phase in which the ceramic granules surrounded by inflammatory cells (1Y2 wk), there were increasing signs of vascularization. Osteoneogenesis and – at the same time - resorption of the HA ceramic were observed after the third week. No major difference in comparison to the bovine bone material could be found.

Discussion: Our results confirm the favorable qualities of the new ceramic reported in association with current maxillofacial literature. Conventional HA granules used for mastoid obliteration to date often showed problems with prolonged inflammatory reactions and, finally, extrusions. In contrast to those ceramics, the new material seems to induce more osteoneogenesis and undergoes early resorption probably due to its high porosity. Overall, it is similar to the bovine bone substance tested on the opposite ear in each animal. Further clinical studies may reveal whether NanoBone® can be an adequate material for obliterating open mastoid cavities in patients.

Harms C, Helms K, Taschner T, Stratos I, Gerber T, Lenz S, Vollmar B, Mittlmeier T

Histomorphometric and micro-CT analysis of the osteoneogenic capacity in the metaphysis of the sheep after implantation of nanocrystalline bone grafting substitute NanoBone™

Chirurgisches Forum 2008, Band 37; 253:255

Abstract: Autologous cancellous bone transplantation is today the gold standard to substitute large bone defects. However a high rate of transplant morbidity and a limited transplant availability reduce the use of this therapy. The synthetic material NanoBone™ (hydroxylapatite nanocrystallines embedded in a porous silica gel matrix) has been shown to have in vivo osteoconductive properties on desmal bone. Up to now positive clinical experience with NanoBone™ has been reported after ventral body fusions and in mandibular surgery. Goal of our study was to examine in vivo the applicability of NanoBone™ on long bones, using a standardized bone defect model for the sheep tibial metaphysis. Therefore we used 18 full-grown sheep and milled a 1.05 cm³ standardized defect under the articular surface of the medial tibia condyles on both hind legs. The defect on the right was filled up with NanoBone™ and the defect on the contralateral leg was left empty. Because of the compact cancellous structure of the tibia head and the integrity of the lateral condyles, a further stabilization was not needed. The defect was partially loaded under full load of the hind leg. Animals were sacrificed after 6 weeks (n = 6), 12 weeks (n = 6) and 26 weeks (n = 6). Operated bones were explanted, taken into account their periosteal integrity. Specimens were macroscopically and microscopically analyzed. Radiographic analysis was performed by means of X-ray, macro- and micro-CT. By compiling the volume distribution and radiographic density from the micro-CT data, a minimal value for osteoneogenesis was calculated. Upon decalcification samples were histomorphometrically analyzed.

The histological and radiological analysis of the defect on the left control-side showed no bone formation after 6, 12 and 26 weeks. In contrast, the micro-CT analysis of the right with NanoBone™ filled defect showed a 55 % volume fraction of structures with bone density. Furthermore the quantitative histochemical analysis of 6 weeks revealed an osteoneogenesis of 22 % and of 12 weeks 34 %. After 12 weeks the micro-CT analysis showed an increase of the structures with bone density to 72 % and after 26 weeks to 74 %. HE-sections demonstrated multinucleated giant cells on the surface of the biomaterial and resorption lacunae, indicating resorption by osteoclasts. In conclusion NanoBone™ is a highly potent bone replacement material with osteoconductive properties in sheep, supporting the potential use of NanoBone™ also in humans.

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Schrodi I, Abshagen K, Gerber T, Vollmar B

In vivo analysis of biocompatibility and vascularization of the synthetic bone grafting substitute NanoBone™

Chirurgisches Forum 2008, Band 37; 251:252

Abstract: One of the major challenges in tissue engineering of bone substitutes are adequate vascularization and biocompatibility of the implant. Thus, the temporal course of neovascularization and the microvascular inflamma-

tory response of implants of **NanoBone™** (fully synthetic nanocrystalline bone grafting material) were studied in vivo by using the dorsal skinfold chamber model. Angiogenesis, microhemodynamics and leukocyte-endothelial cell interaction were analyzed repetitively after implantation in the center and in the border zone of the implant up to 15 days. Both **NanoBone™** granules and plates exhibited high biocompatibility comparable to that of cancellous bone, as indicated by a lack of venular leukocyte activation after implantation. In both synthetic **NanoBone™** groups, signs of angiogenesis could be observed even at day 5 after implantation whereas granules showed higher functional vessel density compared with **NanoBone™** plates. The angiogenic response of the cancellous bone was markedly accelerated in the center of the implant tissue. Histologically, implant tissue showed an ingrowth of vascularized fibrous tissue into the material combined with an increased number of foreign-body giant cells. In conclusion, **NanoBone™**, particularly in granular form, shows high biocompatibility and high angiogenic response, thus improving the healing of bone defects. Our results underline, that beside the composition and nanostructure, also the macro-structure is of importance for the incorporation of the biomaterial by the host tissue.

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Human periodontal fibroblast response to a nanostructured hydroxyapatite bone replacement graft in vitro

Archives of Oral Biology 2008; 53:683-689

Objective: The efficacy of nanostructured hydroxyapatite (NHA) for the treatment of osseous defects has been demonstrated in recent studies, even though the underlying biological mechanism is still poorly known. This study examined the alterations in cellular adhesion and mitogenic responses in human periodontal ligament (PDL) cells treated with a novel nanostructured hydroxyapatite bone graft substitute and characterized associated changes in cellular signalling pathways.

Methods: Cultured PDL cells were stimulated with NHA in a surface coated form. Proliferation was determined by bromodeoxyuridine (BrdU) incorporation and cell adhesion was analysed by a colorimetric assay. In order to understand altered adhesion properties of PDL fibroblasts their integrin profile was analysed and the phosphorylation status of focal adhesion kinase (FAK) and $\alpha 5 \beta 1$ integrin was determined by immunoblotting. In order to understand the signalling mechanisms of increased cell proliferation of PDL cells caused by NHA, the phosphorylation status of the serine/threonine protein kinase Akt, of the signal regulated kinases ERK1/2 and of the epidermal growth factor receptor (EGFR) was analysed by western blot using phospho-specific antibodies.

Results: The results indicated that NHA is a strong stimulator of PDL cell attachment and proliferation. Mechanistically, $\alpha 5 \beta 1$ integrin-mediated cellular adhesion of PDL fibroblasts, which resulted in altered phosphorylation and activation levels of FAK. Proliferation mediated by NHA was mechanistically caused by activation of the epidermal growth factor receptor (EGFR) pathway and its downstream targets ERK1/2 and Akt.

Conclusions: In sum, our findings present evidence that $\alpha 5 \beta 1$ integrin-mediated cellular adhesion of NHA to PDL fibroblasts, whereas proliferation was caused by activation of the epidermal growth factor receptor (EGFR) and the MAP kinase (ERK1/2) and Akt pathways.

Stübinger S, Ghanaati SM, Orth C, Booms P, Kirkpatrick C, Sader R

A new nano-structured and synthetic biomaterial promotes reconstruction of alveolar ridge defects after dental trauma: A preliminary report of clinical and animal studies

IADT 2008, Poster

Objectives: The following study was undertaken to clinically evaluate the properties of NanoBone™ (a new entirely synthetic and nano-structured hydroxylapatite based biomaterial) as a grafting material for guided bone regeneration after dental trauma. This study was triggered by initial in vivo analysis of the host-biomaterial-interaction in the subcutaneous implantation model in Wistar-rats.

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Summary: Based on experiences made during the last two and a half years we present data on clinical and radiological findings as well as histologic specimens after augmentative treatment using a nanostructured bone substitute (NanoBone™) in comparison to other bone substitutes.

Comparing reports on quantitative evaluations 9 to 12 months after sinus floor elevation with bovine bone matrix (BBM) or β -TCPs in literature to those specimens where we have used NanoBone™ it becomes evident that NanoBone™-cases show a similar or even higher amount of de novo bone formation after only 2 to 3 months. Correlating those results with time and the impressive contrast in the histological sections which show faster and more complete de novo bone formation when NanoBone™ was used our proposal to a substantial shortcut of therapy protocols gets support.

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Summary: This study presents data found in 14 individuals where bilateral sinus floor elevation had to be performed prior to implant therapy. To evaluate the influence of added autogenous bone on the de novo bone formation in sinuses at randomly chosen sides NanoBone™ was used as bone substitute either without or in a mixture with autogenous bone chips. At the second stage procedure 8 to 14 weeks later bone cylinders were collected for histological and histomorphometric analysis. NanoBone™ without additional bone revealed an average of 39.5 % de novo bone while in the group with both bone and NanoBone™ 40.7 % of new bone could be measured. The difference is of no statistical significance

Meier J, Heine M, Wolf E

Shortening Therapy Protocols by using the Nanocrystalline Bone Substitute NanoBone™ for Sinus Floor Elevations and Augmentation of other Bone Defects

EAO 2007, Barcelona, Poster

Introduction: The preparation of a sufficient bony layer prior or simultaneously to implantations is mandatory for good long term results. Depending on the amount of bone missing this can be achieved by one- or two-staged approaches. Those bone substitutes that were introduced previously are characterized by long healing periods of (6-) 9 to 12 months that are required to gain sufficient bone regeneration (bovine bone matrix as well as β -TCP). The use of nanocrystals leads to a faster turnover due to the enlarged surface. The most important question now

is whether those bone substitutes can guarantee the growth of de novo bone and the remodelling that is necessary for the primary stability of endosseous implants. We performed several studies to prove the bony integration of NanoBone™ and estimate the time required for the generation of bone that provides us with the reliable conditions to insert implants at shorter therapy intervals compared to other bone substitutes.

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Presence of osteoclast-like cells in the subcutaneous tissue of Wistar rats: in vivo Biocompatibility analysis of a synthetic HA and SiO₂ matrix

21st European Conference of Biomaterials, Poster (Brighton, UK, 9-13th September 2007)

Sol-gel technology results in a variety of biomaterial surfaces and leads to an enlargement of the interface between the biomaterial and the peri-implant tissue. The bone substitute NanoBone® consists of nanocrystalline hydroxyapatite embedded in a highly porous matrix of silica gel and is produced in a sol-gel process at a temperature of < 700°C.

Evaporation leads to the formation of small pores (Ø 5-100 µm). The crystallites are loosely packed and held together by SiO₂ which connect the HA crystals and leads to nano-pores (Ø10-20 nm).

This interconnective porosity is the characteristic of this new biomaterial and assumed to be responsible for induction of the new bone. Up to now there is no in vivo investigation analysing the biodegradation of this biomaterial and its influences on peri-implant cells in the subcutaneous implantation model. Using histological and histochemical methods our aim was to analyse the biocompatibility of NanoBone™ and to identify the cells involved in this degradation.

Hebecker R, Sola S, Mann S, Buchholz K, Piek J

Lumbar Interbody Fusion with a New Nanostructured HA Bone Substitute (NanoBone™) – A Prospective Clinical and CT Study with 15 Patients

Biospine 2, 2nd International Congress Biotechnologies for Spinal Surgery, Poster (Leipzig, Germany, September 20th-22nd, 2007)

In spinal surgery limited availability of autologous bone graft and donor site morbidity are challenge for bone substitutes becoming more and more important. A variety of materials have been introduced for intervertebral cage filling. Therefore bioceramics as well as osteogenetic growthfactors have mainly been used in recent studies. As a new hydroxyapatite (HA)-based bone substitute (NanoBone™) has already been established for use in craniomaxillofacial surgery with promising results we initiated a prospective study to prove its sufficient bone graft potential for interbody fusion.

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Meier J, Wolf E

Histomorphological and immunohistological findings after sinuslift procedures

Osteology Symposium (Monaco – May 10th–12th, 2007), Poster

The production of bone substitute with nanostructure has significantly enhanced the bone (re)modelling by better absorption of osteogenetic substances at the enormously surface and the high porosity. The nanostructured new bone substitute NanoBone™ provides a material that can be used for sinus floor elevations succeeded by the insertion of dental implants after a healing period of about 3 months in two stage procedures or the beginning of functional loading in one stage procedures after only 3 to 4 months.

Here the cellular ingrowth and formation of new bone surrounding the particle of this bone substitute is presented in histologic sections and compared to others with special regard to the rate of bone formation which was evaluated not only by morphologic appearance but by histomorphometry.

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(Osteoinductivity of a synthetic bone replacement material)

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A New Osteoinductive Bone Replacement Material

International Proceedings, XVIII Congress of the European Association for Cranio-Maxillo facial Surgery, Barcelona (Spain), September 12-15, 2006, 19-22

Introduction Hydroxyapatite (HA) is the main component of bone and an important material used for bone substitutes. Conventional HA ceramics are osteoconductive, but poorly degradable. A new HA-silica-matrix is presented being highly osteoconductive and fully biodegradable. Ectopic bone formation was induced when implanted subcutaneously into fatty tissue proving osteoinductive properties of the new biomaterial.

Henkel KO, Gerber T, Lenz S, Gundlach KH, Bienengräber V

Macroscopical, histological, and morphometric studies of porous bone-replacement materials in minipigs 8 months after implantation

Oral Surg Oral med Oral Pathol Oral Radiol Endod 2006; 102:606-13

Objective: The aim of this investigation was to test the induction of bone formation and biodegradation of different biomaterials based on calcium phosphate (CaP). Up to now, hydroxyapatite and β -tricalcium phosphate ceramics have routinely been sintered at temperatures of 1300°C. The new CaP biomaterials tested are fabricated by a sol-gel process at only 700°C.

Study design: Critical-size defects (>5 cm³) in the mandible of 15 adult Goettingen minipigs were filled with 1 of the 2 new types of CaP biomaterials, or with 1 of 2 well-known old-type ceramics, or with a gelatin sponge (in the control group). Macroscopical, histological, and morphometric examination of the former defect areas were made 8 months postoperatively.

Results: Eight months after implantation of the new CaP biomaterials, complete bone formation was observed in the defect area, and at the same time, the foreign material was resorbed almost completely. After implantation of the classical types of ceramics, only incomplete bone formation and a lesser resorption rate of the foreign bodies were noted. The difference in the bone formation rate was significant: more than 93 % for the new CaP biomaterials versus less than 58 % for the classical types of ceramics ($P < 0.01$).

Conclusion: The biological behavior of the new CaP biomaterials was better than that of the old-type sintered ceramic bone-grafting materials. These new CaP matrices are suitable for filling bone defects and are of interest for dentists, including implantologists, craniomaxillofacial and orthopedic surgeons, as well as traumatologists.

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A new synthetic bone replacement material with osteoinductive properties – in vivo investigations

BIOmaterialien 7 (S1),2006;80

Introduction: Hydroxyapatite (HA) being the main component of bone is an important material used for bone substitutes. Conventional HA ceramics are osteoinductive, but poorly degradable. A new HA matrix is presented being highly osteoconductive and at the same time fully biodegradable. Ectopic bone formation was induced when implanted subcutaneously into fatty tissue proving the osteoinductive properties of the new material.

Gerber T, Holzhüter G, Götz W, Bienengräber V, Henkel KO, Rumpel E
Nanostructuring of Biomaterials – A Pathway to Bone Grafting Substitute
Eur J Trauma 2006;32:132-40

Background The bone substitute **NanoBone™** consists of nanocrystalline hydroxyapatite embedded in a highly porous matrix of silica gel. It promotes the healing of bone defects and is degraded by osteoclasts during bone remodeling. The present study investigates the interactions of **NanoBone™** with bone tissue.

Methods: Granules of **NanoBone™** were implanted in defects of critical size in the mandible of minipigs. Samples were taken after 5 and 10 weeks and demineralized. The composition of the implanted granules was analyzed by means of transmission and scanning electron microscopy and EDX. Enzyme and immunohistochemistry was used to investigate organic components of **NanoBone™** granules that arised after implantation in the host.

Results: EDX demonstrated that 5 weeks after implantation the silica gel was degraded and replaced by an organic matrix. Ultrastructurally, the matrix appeared amorphous with only single collagen fibrillae. PAS-staining indicated the presence of carbohydrates. Immunohistochemically, the bone proteins osteopontin, osteocalcin and BMP-2 were found as constituents of the new matrix. Alkalic phosphatase activity was located in osteoblasts and newly formed bone on **NanoBone™** and focally in particles. Osteoclasts with ruffled borders, sealing zones, and acid phosphatase activity were situated in resorption lacunae at granule surfaces not covered by new bone.

Conclusions: In vivo, the silica gel of **NanoBone™** is replaced by bone matrix glycoproteins with known functions in attraction, adhesion, and differentiation of bone cells as osteoblasts and osteoclasts. We assume that the deposition of these molecules supports the early phase of **NanoBone™** degradation by osteoclasts and promotes the production of new bone tissue.

Dietze S, Bayerlein T, Proff P, Hoffmann A, Gedrange T
The ultrastructure and processing properties of Straumann Bone Ceramic and NanoBone™.
Folia Morphol (Warsz). 2006 Feb;65(1):63-5.

The ultrastructure, fundamental chemistry, and processing modes of fully synthetic bone grafting materials are relevant to the reconstruction of osseous defects. Rapid progress in the profitable market of biomaterials has led to the development of various bone substitutes. Despite all these efforts, an ideal and full substitute of autologous bone is not yet in sight. With regard to anorganic calcium phosphate ceramics, Straumann Bone Ceramic and **NanoBone™** are compared. These have a similar composition and are osteoconductive, which indispensably requires contact with well-vascularised bone.

Gerike W, Bienengräber V, Henkel KO, Bayerlein T, Proff P, Gedrange T, Gerber T
The manufacture of synthetic non-sintered and degradable bone grafting substitutes.
Folia Morphol (Warsz). 2006 Feb;65(1):54-5.

A new synthetic bone grafting substitute (**NanoBone™**, ARTOSS GmbH, Germany) is presented. This is produced by a new technique, the sol-gel-method. This bone grafting substitute consists of nanocrystalline hydroxyapatite (HA) and nanostructured silica (SiO₂). By achieving a highly porous structure good osteoconductivity can be seen. In addition, the material will be completely biodegraded and new own bone is formed. It has been demonstrated that **NanoBone™** is biodegraded by osteoclasts in a manner comparable to the natural bone remodelling process.

Rumpel E, Wolf E, Kauschke E, Bienengräber V, Bayerlein T, Gedrange T, Proff P.

The biodegradation of hydroxyapatite bone graft substitutes in vivo.

Folia Morphol (Warsz). 2006 Feb;65(1):43-8.

Hydroxyapatite (HA) ceramics are widely used for bone reconstruction. They are osteoconductive and serve as structural scaffolds for the deposition of new bone. Generally, scaffold materials should be degradable as they affect the mechanical properties of the reconstructed bone negatively. Degradation by osteoclasts during the bone remodelling process is desirable but often does not take place. In the current study we analysed by light microscopy the degradation of two granular HA implants in critically sized defects in the mandibula of Goettingen mini-pigs five weeks after implantation. Bio-Oss consists of sintered bovine bone and NanoBone™ is a synthetic HA produced in a sol-gel process in the presence of SiO₂. We found that both biomaterials were degraded by osteoclasts with ruffled borders and acid phosphatase activity. The osteoclasts created resorption lacunae and resorptive trails and contained mineral particles. Frequently, resorption surfaces were in direct contact with bone formative surfaces on one granule. Granules, especially of NanoBone™, were also covered by osteoclasts if located in vascularised connective tissue distant from bone tissue. However, this usually occurred without the creation of resorption lacunae. The former defect margins consisted of newly formed bone often without remnants of bone substitutes. Our results show that the degradation of both biomaterials corresponds to the natural bone degradation processes and suggest the possibility of complete resorption during bone remodelling.

Kauschke E, Rumpel E, Fanghänel J, Bayerlein T, Gedrange T, Proff P.

The in vitro viability and growth of fibroblasts cultured in the presence of different bone grafting materials (NanoBone™ and Straumann Bone Ceramic).

Folia Morphol (Warsz). 2006 Feb;65(1):37-42.

Different clinical applications, including dentistry, are making increasing demands on bone grafting material. In the present study we have analysed the viability, proliferation and growth characteristics of fibroblasts cultured in vitro together with two different bone grafting materials, NanoBone™ and Straumann Bone Ceramic, over a period of 24 and 28 days respectively. Viability was measured at least every 72 hours by using the alamarBlue assay, a test that measures quantitatively cell proliferation and viability but does not require cell fixation or extraction. After one week of culture fibroblast viability was as high as in controls for both grafting materials and remained high (> 90 %) for the duration of the experiment. Cell growth was evaluated microscopically. Scanning electron microscopy revealed a dense fibroblast growth at the surface of both bone grafting materials after three weeks of in vitro culture. Generally, our in vitro analyses contribute to further insights into cell - scaffold interactions.

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Chuchracky N

NanoBone® Augmentation Material and Bego Semados® - S-Implants: A Powerful Combination for Today's Dental Implantology Applications?

implants 1_2006:06-09

NanoBone® Bone Augmentation Material

The special properties of the all-new NanoBone™ Augmentation material derive from the material's nanostructure, which means that identical chemical compounds can have completely different properties. NanoBone™ has a very large surface area relative to its volume. Numerous studies have shown that the surface properties of bone augmentation materials are of decisive importance for the development of optimal biological activity and the formation of new bone (osteogenesis). The deposition of hydroxyapatite in a SiO₂ structure forms the basis of the material. Silicon dioxide molecules are used to achieve the additional effect of "fixing" proteins in the surface. Silicon dioxide is particularly important here because it stimulates the formation of collagen and bone.

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Henkel KO, Gerber Th, Dörfling P, Gundlach KH, Bienengräber V

Repair of bone defects by applying biomatrices with and without autologous osteoblasts

Journal of Cranio-maxillofacial Surgery (2005) 33, 45-49

Question: Is it possible to stimulate osteoconduction and osteogenesis to improve bone formation in critical-size defects in order to avoid bone grafting? Material and methods: Full thickness, critical-sized defects were created in the anterior mandible of 16 adult mini-pigs. The defects were filled with a new bioactive matrix (60 % hydroxyapatite and 40 % β -tricalciumphosphate), produced by an innovative low temperature sol-gel process (120 °C). The biomatrix was tested alone and in combination with cultured autologous osteoblasts. In a control group, periosteum was the only bone producing source. Five weeks postoperatively, the animals were sacrificed and the defects analysed macroscopically, histologically and radiographically.

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Comparison of a new kind of calcium phosphate formula versus conventional calciumphosphate matrices in treating bone defects – A long-term investigation in pigs

Key Engineering Materials Vols. 284-286 (2005) pp. 885-888

In clinical practice arises an increasing need for bone substitute materials. The main inorganic part of bone is the hydroxyapatite (HA). A new hydroxyapatite formula was created by a sol-gel-process at low temperature level [4]. The aim of this investigation was to test the biodegradation and the induction of bone formation by this new material and to compare these versus conventional fabricated HA and β -TCP. 30 one-year-old Goettingen minipigs were divided into five groups. Critical size defect ($>5 \text{ cm}^3$) in the mandible was treated differently in all 5 groups: group I- filling with pure HA, which was fabricated by sol-gel-technique, group II- control, only gelatinous material was given, group III- conventional β -TCP [Cerasorb®], in group IV- conventional HA [Endobone®] and in group V [Targobone®], a non denatured bovine collagen matrix was used. Macroscopical and microscopical investigations of the former defects were made eight months postoperatively. The bone formation was superior in the sol-gel-HA-group (group I) in comparison with the control groups (group II) and the conventional fabricated ceramics groups (III and IV). In the sol-gel-HA group, the biodegradation of this new biomaterial was considered to very good with a resorption rate of more than 98 %; eight months postoperatively. In this group complete bone formation was seen in former defects. In the control group, only an incomplete bone formation with 48.4 % of the defect area was noted. This difference was significant ($p < 0,001$).

A less bone formation was also observed in group III and IV with 57.6 % and 56.9 %. The bovine non-denatured collagen matrix (group V) leads to only 20 % of new formed bone. The new calcium phosphate formula made by a sol-gel method seems to be superior and suitable for filling bone defects.

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Background: Up to now hydroxyapatite (HA) and β -tricalciumphosphate (β -TCP) ceramics have been routinely sintered at temperatures between 1100° and 1500°C. Our new calcium ceramic is fabricated by a sol-gel process at 200°C. The aim of this investigation was to test the biodegradation of and the induction of bone formation by this material.

Material and methods: Eighteen 1-year-old Goettingen minipigs were divided into three groups. Critical size defects ($>5 \text{ cm}^3$) in the mandible were treated differently in all three animals (group 1: filling with 40 % β -TCP plus 60 % HA, group 2: pure HA was applied, group 3 served as controls: only gelatinous material was given). Macroscopic and microscopic investigations of the former defects were made 8 months postoperatively.

Results: In groups 1 and 2 biodegradation of more than 93 % of the new calcium phosphate formula was found 8 months postoperatively and considered to be very good. No difference was observed between pure HA (group 2) and the combination of HA and β -TCP (group 1). In both groups complete bone formation was seen macroscopically in the former defects. In the control group only incomplete bone formation with 48.4 % of the defect area was noted. This difference was significant ($p<0.001$).

Discussion: The new calcium phosphate formula made by a sol-gel method at 120°C seems to be suitable for filling bone defects and is of interest for orthopedic surgery, traumatology, craniomaxillofacial surgery, and dentistry.

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