WHAT'S NEW DEUTSCH

exoplan
3.1 Rijeka

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exocad

Tooth: 36 (FDI) | Patient: Doe, John | Project: 2020-12-01_00001-001 (Lower jaw)

4.6 mm Master Cylinder H6

Article Number	CGS-QMC-10
Inner Diameter Ø	5.11 mm
Height H	5.99 mm
To platform: 10.53 mm	
Distance from sleeve top to bone level: 10.53 mm	

Biocorizons® - Guided Surgery Kit CGS4000

Step	Drill Size	Length
1	Ø 4.6 mm	L 24 mm
2	Ø 2 mm	L 24 mm
3	Ø 2.5 mm	L 24 mm
4	Ø 3.2 mm	L 24 mm
5	Ø 3.7 mm	L 24 mm
6	Ø 4.1 mm	L 24 mm
7	Stop Position SP3	

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
Implant Positioning | 32 Search anything... | **Guided** | **Favorites**

Implant systems (83)

Length (L)	Ø 3.5 mm RB	Ø 3.75 mm RB	Ø 4.0 mm RB	Ø 4.5 mm RB	Ø 5.0 mm WB
6 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
8 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
10 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
12 mm	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
14 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
16 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
18 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

32 straumanngroup

BLX - Bone Level X Roxolid® SLActive®

• Ø 3.75 mm RB, 12 mm

Select surgical component

Select prosthetic component

42 Tooth 42

Select implant

Manufacturer: Straumann® Group

System: BLX - Bone Level X Roxolid® SLActive®

Implant Analog: BLX, Ø 3.75 mm RB, SLActive® 12 mm, Roxolid®

Article No. 061.4312

Length: 12 mm | **Platform Ø: 3.5 mm** | **Body Ø: 3.75 mm**

Comment: Material: Roxolid®, Surface: SLActive®

← BACK | NEXT → | You can now place the implant with the mouse. | Sleeves / Kits > | Prosthetics > | Confirm

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Loaded series: F, Series date: 5/15/2021

Study: Dual-Scan-Protocol, Example

Series: Patient with markers in both jaws (407 Images)

Modality: CT
Description: Patient with markers in both jaws
No. of images: 407
Image type: AXIAL
Series date: 5/15/2021

Buttons: Select file set, Load series, Discard series

Approval of Planning

The implant planning is now completed. You are now kindly asked to approve your planning and confirm the items below before the implant planning report and data files are generated.

Patient Information		
Name	Project	DICOM
[Blurred]	[Blurred]	[Blurred]
Birth date		


- I agree to study the implant planning report before commencing surgery and verify the provided information with regard to the surgery.
- I confirm that I have performed the planning with due care and that I am satisfied with the medical and clinical results of the planning. I chose the most appropriate implant solution and observed safety margins to any existing dentition, existing and intended dental restorations, and anatomical structures.

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Buttons: Add all images to report, Save all images, Delete all images

Buttons: Capture main view only, Capture manually, Capture all views, Load from file

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