

DryDent®

By DIRECTA

Directa's products are developed, tested and evaluated by our selected team of highly-qualified and renowned Swedish and international dental professionals with the aim of making life easier for the dental professionals in the daily clinical work.

Design By Dentists™

Directa's products are developed, tested and evaluated by our selected team of highly-qualified and renowned Swedish and international dental professionals with the aim of making life easier for the dental professionals in the daily clinical work.

DIRECTA
directadentalgroup.com



DIRECTA AB
P.O. Box 723,
194 27 Upplands Väsby, SWEDEN
+46 8 506 505 75

DIRECTA Inc
64 Barnabas Road, Unit 3
Newtown, CT 06740, United States
203.491.2273

DIRECTA GERMANY
Hammergutstr. 11
92245 Kümmersbruck, GERMANY

www.directadental.com
info@directadental.com

2016-02

DirectaDentalGroup

DIRECTA **topdental** **ORSING**
directadentalgroup.com directadentalgroup.com directadentalgroup.com

DIRECTA
directadentalgroup.com

4 reasons to use DryDent® Parotid

#1

Super Absorbent

DryDent Parotid is placed directly over the open ducts of the parotid gland, maximizing saliva absorption throughout the procedure.

#2

Dry Dentistry

DryDent Parotid's capacity to absorb unwanted moisture quickly provides an optimal dry treatment area every time.



#3

Soft and Gentle

DryDent Parotid's new design with soft edges is made of smooth materials that provide maximum patient comfort.

#4

Improved Properties

DryDent Parotid's unique new foam material without binders maintains its form, keeping the absorbed saliva contained even when fully saturated.

DryDent® Parotid. For optimal moisture control.

Absorbing and Protecting!

The new DryDent Parotid introduces new features for controlling saliva produced by the parotid glands.

DryDent Parotid is placed directly over the open ducts of the parotid gland where saliva is optimally absorbed, allowing for dryer working areas. Its super absorbent foam material traps unwanted moisture, facilitating dryer fields for impressiontaking, cementing and filling procedures.

DryDent Parotid's super absorbent capacity keeps moisture control even when it is fully saturated and it is made of smooth soft materials that provide maximum patient comfort throughout any treatment.

Available in Large and Small packages of 50 pieces.

4 reasons to use DryDent® Sublingual

#1 The First Sublingual Super Absorbent

DryDent Sublingual is the only super absorbent, designed to collect unwanted moisture under the tongue.

#2 Dry Dentistry

DryDent Sublingual absorbing properties collect saliva quickly, providing an optimal dry area for a successful treatment every time.



#3 Soft and Gentle

DryDent Sublingual reduces the swallowing reflexes due to its smooth comfortable design which provides maximum patient comfort.

#4 Improved Comfort

DryDent Sublingual improves patient comfort when noisy and painful saliva ejectors can be placed over the DryDent rather than over soft tissue.



DryDent® Sublingual. Unique saliva control.

Super Absorbents have reached a new level with the introduction of DryDent Sublingual.

A unique new product which collects and controls the saliva produced by the sublingual glands as well as the submandibular glands.

It is placed under the tongue for the optimal collection of unwanted moisture under the sublingual area. Its super absorbent qualities provide a super dry oral environment during a variety of procedures.

Patient comfort has been improved with its smooth soft material design that reduces uncomfortable swallowing reflexes. Saliva ejectors can be placed over the DryDent Sublingual, minimizing uncomfortable suction of soft tissue. DryDent Sublingual's super absorbent capacity keeps moisture control even when it is fully saturated.

Available in Large package of 40 pieces and Small packages of 50 pieces.

Product Numbers

DryDent® Parotid

Prod No	Product Description	Quantity
416073	DryDent® Parotid Small Size: 40x36x2 mm Expand Size: 40x36x10 mm	50 pcs/package
416072	DryDent® Parotid Large Size: 50x43x2 mm Expand Size: 50x43x10 mm	50 pcs/package

DryDent® Sublingual

Prod No	Product Description	Quantity
416069	DryDent® Sublingual Small Size: 30x50x2 mm Expanded Size: 30x50x10 mm	50 pcs/package
416068	DryDent® Sublingual Large Size: 38x60x2 mm Expanded Size: 38x60x10 mm	40 pcs/package

