

Individuelle Abutments gestern vs. heute

Dr. Ali-Reza Kebai

Implantologie Journal 2/2011

Literatur

1. Adams MW. Computer-designed and milled patient-specific implant abutments. *Dent Today* 2005; 24(6): 80-3
2. Ganz S. Computer-milled patient-specific abutments: incredibly quality with unprecedented simplicity. *Implantology* 2003 2003: 37-44
3. Ganz S. Finally, a 'win-win' solution: Increasing accuracy while saving time, money with computer-milled abutments. *Dental Economics* 2005; May
4. Garg A. The Atlantis Components abutment: simplifying the tooth implant procedure. *Dental Implantology Update* 2002;September
5. Kerstein RB, Osorio J. Utilizing computer-generated duplicate titanium custom abutments to facilitate intraoral and laboratory implant prosthesis fabrication. *Pract Proced Aesthet Dent* 2003; 15(4): 311-4
6. Nazarian A. Easier implant restoration: CAD/CAM generated implant abutments. *Contemporary Esthetics* 2007; February: 44-48
7. Schneider A, Kurtzman GM. Computerized milled solid implant abutments utilized at second stage surgery. *Gen Dent* 2001; 49(4):416-20
8. Whitesides L. Evaluation of the Atlantis abutment in implant restoration. *Inside Dentistry* 2006; September: 98-99
9. Watkin A, Kerstein RB. Improving darkened anterior peri-implant tissue color with zirconia custom implant abutments. *Compend Contin Educ Dent* 2008; 29(4): 238-40, 42