

Ausgabe: ZWP Spezial 4/16

Thema: CED-IADR 2015: Positive In-vitro-Daten zu EQUIA Forte

Autoren: Dr. Falk Schwendicke

Literatur:

1. Basso M et al. (2014): 48-Months, Multicentre, Clinical Evaluation on 304 Glass Ionomer Permanent Restorations. Abstract 192686, IADR-PER Dubrovnik 2014.
2. Diem V et al. (2014): The effect of a nano-filled resin coating on the 3-year clinical performance of a conventional high-viscosity glass-ionomer cement. Clin Oral Investig 18(3):753-759.
3. Friedl K et al. (2011): Clinical performance of a new glass ionomer based restoration system: A retrospective cohort study. Dent Mater 27(10):1031-1037.
4. Gurgan S et al. (2014): 60-Month Clinical Performance Of A Glass-Ionomer Restorative System. Abstract 89, IADR Cape Town 2014.
5. Klinke T et al. (2013): Randomized Clinical Trial in the field: Longevity after 24 months. J Dent Res 92 (Spec Iss B), Abstract 3.
6. Klinke T et al. (2015): 48 months clinical performance of two current glass-ionomer systems in a field study. A RCT in field. ConsEuro 2015.
7. GC Corporation: Research & Development.
8. Schwendicke F et al. (2015): Pulpal remineralisation of artificial residual caries lesions in vitro. Caries Research 49:591-594.