



Dosieranleitung

	DesNet+ TR-3 BIB forte eco	DesNet+ Alpron	DesNET+ AlproJet-DD BIB forte eco	BIB forte eco AlproSol	BIB forte eco AlproSol	DesNET+ AlproSol AlproJet-D AlproJet-W TarClean
Gebrauchslösung	0,5%	1,0%	2,0%	3,0%	4,0%	5,0%
1 Liter	5 ml	10 ml	20 ml	30 ml	40 ml	50 ml
2 Liter	10 ml	20 ml	40 ml	60 ml	80 ml	100 ml
3 Liter	15 ml	30 ml	60 ml	90 ml	120 ml	150 ml
4 Liter	20 ml	40 ml	80 ml	120 ml	160 ml	200 ml
5 Liter	25 ml	50 ml	100 ml	150 ml	200 ml	250 ml
6 Liter	30 ml	60 ml	120 ml	180 ml	240 ml	300 ml
7 Liter	35 ml	70 ml	140 ml	210 ml	280 ml	350 ml
8 Liter	40 ml	80 ml	160 ml	240 ml	320 ml	400 ml

REF 803501.18

