

**Ausgabe:** ZWP Zahnarzt Wirtschaft Praxis 12/11, S. 90-94

**Thema:** Vollkeramische CAD/CAM Materialien

**Autor:** Rusell A. Giordano, DMD, DMSc.

### **Literatur**

1. Seghi RR, Sorensen JA: Relative flexural strength of six new ceramic materials. *Int J Prosth* 8:239-246, 1995.
2. Giordano R, Kanchanatawewat K, Asvanund P, Nathanson D: Flexural strength evaluation of ceramics for Celay restorations. *J Dent Res Special Issue* 75:125, 1996.
3. Krejci I, Lutz F, Reimer M: "Wear of CAD/CAM ceramic inlays: restorations, opposing cusps, and luting cements." *Quintessence Int* 1994 Mar;25(3):199-207.
4. Abozenada, B, Pober, R, Giordano, R "In-vitro wear of restorative dental materials" Abstract# 1693, *J. Dent. Res., Special Issue, VOL 81, MAR 2002*.
5. McLaren, E., Abuzenada, B, Pober, R, Giordano, R; "Material Testing and Layering Techniques of a New Two-Phase All-Glass Veneering Porcelain for Bonded Porcelain and High-Alumina Frameworks" *QDT 2003 (Vol 26)*.
6. McLaren, E., Giordano, R. " Zirconia-Based Ceramics: Material Properties, Esthetics, and Layering Techniques of a New Veneering Porcelain, VM9, High-Alumina Frameworks" *QDT 2005 (Vol 28)*.
7. Tinschert J, Zwez D, Marx R, Anusavice KJ. Structural reliability of alumina-, feldspar-, leucite-, mica- and zirconia-based ceramics. *J Dent.* 2000 Sep;28(7):529-35.
8. Martin N, Jedyakiewicz NM. Clinical performance of CEREC ceramic inlays: a systematic review. *Dent Mater.* 1999 Jan;15(1):54-61.
9. Mörmann, W.H. (ed): *CAD/CIM in Aesthetic Dentistry, Cerec 10 Year Anniversary Symposium*, Quintessence, Chicago, 1996.
10. Bindl A, Mormann WH.; Survival rate of mono-ceramic and ceramic-core CAD/CAM-generated anterior crowns over 2-5 years., *Eur J Oral Sci.* 2004 Apr;112(2):197-204.
11. Posselt A, Kerschbaum T.; Longevity of 2328 chairside Cerec inlays and onlays., *Int J Comput Dent.* 2003 Jul;6(3):231-48.
12. Studer S, Lehner C, and Scharer P: Seven year results of leucite reinforced glass ceramic inlays and onlays. *J Dent Res Special Issue* 77:802, 1998.
13. Lehner C, Studer S, Brodbeck U, and Scharer P: Seven year results of leucite reinforced glass ceramic crowns. *J Dent Res Special Issue* 77:802, 1998.
14. Fasbinder DJ. Clinical performance of chairside CAD/CAM restorations. *J Am Dent Assoc.* 2006 Sep;137 Suppl:22S-31S.
15. Christensen, G. *CRA News*, April 2006.
16. Lee, Dal Ho " Accuracy of Milled Porcelain, Ceramic and Composite Resin Materials", Masters Thesis Boston university, Goldman School of Dental Medicine, 2000.
17. crowns. A SEM investigation" Abstract#2641, *J. Dent. Res., 2003 Special Issue MAR*.

18. Bindl A et. Al., "Survival of ceramic CAD/CAM crowns bonded to preparations with reduced macroretention geometry." *I J P.* 2005 May- Jun;18(3):219-24
19. Pilo R, Brosh T, Chweidan H: Cusp reinforcement by bonding of amalgam restorations. *J Dent* 26(5-6):467 – 472, 1998.
20. Magne P, Douglas WH: Dentin bonding optimization and biomimetic recovery of the crown. *Int. J. Prosthodont* 12(2):111-121, Mar-Apr 1999.
21. Attia A, Kern M.; Fracture strength of all-ceramic crowns luted using two bonding methods., *J Prosthet Dent.* 2004 Mar;91(3):247-52.