

CEREC und Feldspatkeramik: Höchste Ästhetik in nur einer Sitzung

Dr. Julián Conejo Gutiérrez/Costa Rica

Dental Tribune Austrian/Swiss Edition 9-2015

Literatur

- BINDL A, MÖRMANN WH: Clinical and technical aspects of the CEREC In-Ceram crown. In: Mörmann, WH (ed): CAD/CAM in Aesthetic Dentistry. CEREC 10 Year Anniversary Symposium. Chicago: Quintessence 1996: 441-462.
- BITZ M, MÖRMANN WH: Documentation of CEREC crowns via electronic communication. In: Mörmann, WH (ed): CAD/CAM in Aesthetic Dentistry. CEREC 10 Year Anniversary Symposium. Chicago: Quintessence 1996: 495-516.
- LAMPE K, LÜTHY H, MÖRMANN WH: Fracture load of all-ceramic computer crowns. In: Mörmann, WH (ed): CAD/CAM in Aesthetic Dentistry. CEREC 10 Year Anniversary Symposium. Chicago: Quintessence 1996: 463-482.
- DEVIGUS A: Die Cerec-2-Frontzahnkrone. In: Dental Magazin 1997; 3: 38-41.
- LAMPE K, LÜTHY H, MÖRMANN WH, LUTZ F: Bruchlast vollkeramischer Computerkronen. In: Acta Med Dent Helv 1997; 2: 76-83.
- MÖRMANN WH, BINDL A, LÜTHY H, RATHKE A: Effects of preparation and luting system on all-ceramic computer-generated crowns. In: Int J. Prosthodont 1998; 11: 333-339
- BINDL A, MÖRMANN WH: Clinical evaluation of adhesively placed CEREC Endo-Crowns after 2 years – Preliminary results. In: J Adhesive Dent 1999; 1: 255–265.
- CHEN HY, HICKEL R, SETCOS JC, KUNZLEMANN KH: Effects of surface finish and fatigue testing on the fracture strength of CAD/CAM and pressed ceramic crowns. In: J Prostet Dent 1999; 82: 468-475.
- MALAMENT KA, SOCRANSKY SS: Survival of Dicor glass-ceramic dental restorations over 14 years: Part I. Survival of Dicor complete coverage restorations and effect of internal surface acid etching, tooth position, gender and age. In: J Prosthet Dent 1999; 81: 23-32.
- TOTH RT, MÖRMANN WH: CEREC 2 vollkeramische CAD/CAM Computerkronen – Konstruktion im Front- und Seitenzahnbereich. CAD/CAM Bibliothek Band 1, SFCZ-Verlag Zürich 1999.
- CHAI J, TAKAHASHI Y, SULAIMAN F, CHONG K, LAUTENSCHLAGER EP: Probability of fracture of all-ceramic crowns. In: Int J Prosthodont 2000; 13: 420-424.

- MÖRMANN WH, BINDL A: The CEREC 3 – A quantum leap for computer-aided restorations: Initial clinical results. In: Quintessence Int 2000; 31: 699-712.
- POLANSKY R, ARNETZL G, HAAS M, KEIL C, WIMMER G, LOREN-ZONI M: Restdentinstärke nach 1,2 mm-Stufenpräparation für CEREC-Kronen. In: Int J Comput Dent 2000; 3:243-258.
- BESIMO CE, SPIELMANN H-P, ROHNER H-P: Computer-assisted generation of all-ceramic crowns and fixed partial dentures. In: Int J Comput Dent 2001; 4: 243-262.
- TSOTSOS S: Single-appointment, all-ceramic anterior restorations. In: Int J Comput Dent. 2001; 4: 263-72.
- BINDL A, MÖRMANN WH: An up to 5-year Clinical Evaluation of Posterior In-Ceram CAD/CAM Core Crowns. In: Int J Prosthodont 2002; 15: 451-456.
- MCLAREN EA, TERRY DA: CAD/CAM systems, materials, and clinical guidelines for all-ceramic crowns and fixed partial dentures. In: Compend Contin Educ Dent 2002; 23: 637-641.
- LAWN BR, DENG Y, LLOYD IK, JANAL MN, REKOW D, THOMPSON VP: Materials design of ceramic-based layer structures for crowns. In: J Dent Res 2002; 81: 433-438.
- REICH S, HORNBERGER H: The effect of multicolored machinable ceramics on the esthetics of all-ceramic crowns. J Prosthet Dent 2002; 88: 44-49.
- NAKAMURA T, DEI N, KOJIMA T, WAKABAYASHI K: Marginal and internal fit of CEREC 3 CAD/CAM all-ceramic crowns. In: Int J Prosthodont 2003; 16: 244-248.
- ATTIA A, KERN M: Influence of cyclic loading and luting agents on the fracture load of two all-ceramic crown systems. In: J Prosthet Dent 2004; 92: 551-556.
- OTTO T: Computer-aided direct all-ceramic crowns: preliminary 1-year results of a prospective clinical study. In: Int J Perio Rest Dent 2004; 24: 447-455.
- REICH SM, WICHMANN M, RINNE H, SHORTALL A: Clinical performance of large, all-ceramic CAD/CAM-generated restorations after three years: a pilot study. In: J Am Dent Assoc 2004; 135: 605-612.