

**Ausgabe: ZWP 4/2014**

**Thema: 2-Jahresstudie: Glasionomer-basiertes Restaurationssystem**

**Autoren:** Prof. Dr. Sevil Gurgan, Dr. Esra Firat, Zeynep Bilge Kutuk

## **Literatur**

- <sup>1</sup> Ilie N et al.: Fracture toughness of dental restorative materials. *Clin Oral Investig* 16(2):489–498. (2012)
- <sup>2</sup> Lohbauer U et al.: Strength and wear resistance of a dental glass-ionomer cement with a novel nanofilled resin coating. *Am J Dent* 2011; 24: 124–128. (2011)
- <sup>3</sup> Triana, R et al.: Dentin bond strength of fluoride-releasing materials. *Am J Dent* 7: 252. (1994)
- <sup>4</sup> Forsten L. Fluoride release and uptake by glass-ionomers and related materials and its clinical effect. *Bio-materials* 1998; 19:503-508. (1998)
- <sup>5</sup> Kantovitz K.R. et al.: Inhibition of mineral loss at the enamel/sealant interface of fissures sealed with fluoride- and non-fluoride containing dental materials in vitro. *Acta Odontol Scand* 2006; 64(6):376-383.
- <sup>6</sup> Mjör IA. Glass-ionomer cement restorations and secondary caries: a preliminary report. *Quintessence Int* 1996; 27:171-174. (1996)
- <sup>7</sup> Kielbassa AM, Muller J, Gernhardt CR: Closing the gap between oral hygiene and minimally invasive dentistry: a review on the resin infiltration technique of incipient (proximal) enamel lesions. *Quintessence Int* 2009; 40: 663–681. (2009)
- <sup>8</sup> Lohbauer U et al.: Strength and wear resistance of a dental glass-ionomer cement with a novel nanofilled resin coating. *Am J Dent* 2011; 24: 124–128. (2011)
- <sup>9</sup> Magni E et al.: SEM and microleakage evaluation of the marginal integrity of two types of class V restorations with or without the use of a light-curable coating material and of polishing. *J Dent* 2008; 36: 885–891. (2008)
- <sup>10</sup> Gurgan S et al.: 24-months evaluation of a posterior glass-Ionomer restorative using replication technique. IADR Helsinki, 2012, Abstract 381. Abruf am 12.08.2013 unter: <https://iadr.confex.com/iadr/per12/webprogram/Paper168298.html>
- <sup>11</sup> Gurgan S et al.: 36-months clinical performance of a glass-ionomer restorative system, IADR Seattle, 2013, Abstract 2933. Abruf am 12.08.2013 unter: <http://www.iadr.confex.com/iadr/13iags/webprogram/Paper176766.html>
- <sup>12</sup> Diem VT: The effect of a nano-filled resin coating on the 3-year clinical performance of a conventional high-viscosity glass-ionomer cement. *Clin Oral Investig*. 2013 Jul 7. [Epub ahead of print]