

## Literatur:

1. United Nations Environment Programme (2013). Intergovernmental negotiating committee to prepare a global legally binding instrument on mercury. Fifth session, Geneva, 13–18 January 2013. Abruf am 08.08.2014 unter:  
[http://www.unep.org/chemicalsandwaste/Portals/9/Mercury/Documents/INC5/INC5\\_7\\_asterix\\_final%20report\\_26%2008\\_e.pdf](http://www.unep.org/chemicalsandwaste/Portals/9/Mercury/Documents/INC5/INC5_7_asterix_final%20report_26%2008_e.pdf)
2. Mickenautsch S, Yengopal V (2012). Failure rate of high-viscosity GIC based ART compared to that of conventional amalgam restorations - evidence from a systematic review update. *S Afr Dent J* 2012; 67: 329-31.  
<http://www.ncbi.nlm.nih.gov/pubmed/23951787>
3. Mickenautsch S, Yengopal V, Banerjee A (2010). Atraumatic restorative treatment versus amalgam restoration longevity: a systematic review. *Clin Oral Investig* 2010 Jun;14(3):233-40.
4. Manhart, J, Chen, HY, Hamm, G und Hickel, R (2004): Review of the clinical survival of direct and indirect restorations in posterior teeth of the permanent dentition. *Oper Dent* 29, 481-508
5. Gemeinsame Stellungnahme der Deutschen Gesellschaft für Zahnerhaltung (DGZ) und der Deutschen Gesellschaft für Zahn-, Mund- und Kieferheilkunde (DGZMK) (2005). Direkte Kompositrestaurationen im Seitenzahnbereich – Indikation und Lebensdauer. Abruf am 23.08.2014 unter:  
[http://www.dgzmk.de/uploads/tx\\_szdgzmkdocuments/Kompositrestaurationen\\_im\\_Seitenzahnbereich.pdf](http://www.dgzmk.de/uploads/tx_szdgzmkdocuments/Kompositrestaurationen_im_Seitenzahnbereich.pdf)
6. Diem VTK et al. (2013). The effect of a nano-filled resin coating on the 3-year clinical performance of a conventional high-viscosity glass-ionomer cement. *Clinical Oral Investigations* DOI 10.1007/s00784-013-1026-z
7. Klinke T, Daboul A, Biffar R (2013). EQUIA - RCT in the field: Longevity after 24 months. CED IADR Florenz, 2013, Abstract 3. Abruf am 23.08.2014 unter:  
<https://iadr.confex.com/iadr/ced13/webprogram/Paper179792.html>
8. Gurgan S, Kutuk ZB et al. (2014). 60-Month Clinical Performance Of A Glass-Ionomer Restorative System, IADR General Session, Cape Town, June 25th 2014, Abstract of Oral Session. Abruf am 23.08.2014 unter:  
<https://iadr.confex.com/iadr/14iags/webprogram/Paper190266.html>